

# LarySeal™ Pro

## Laryngeal Mask Airway

Product Data Sheet

**PRODUCT:** LarySeal™ Pro**DESCRIPTION / FUNCTION**

LarySeal™ Pro is a supraglottic airway device (SAD) that quickly facilitates a safe airway and provides gastric access and intubating ability. With a curve that better suits the oropharyngeal anatomy and symmetrical cuff with a pointed distal end, the semi rigid LarySeal™ Pro allows for easy insertion into the patient's airway. Attached to the mask is a cuff inflation line with a color-coded pilot balloon and non-return valve for mask inflation and deflation. The ISO standard 15M connector provides a secure connection to all standard breathing systems.

LarySeal™ Pro can be used as a conduit for intubation using an Endotracheal Tube (ETT). This is further aided by the epiglottic flap (Sizes 3 - 5 only) which helps to lift the epiglottis when pushing through the ETT. The design of the LarySeal™ Pro shell contains a semi rigid diverter inside the space of the ventilatory opening that allows the ETT to be diverted over the suction channel and through the vocal cords.

Additionally, LarySeal™ Pro contains an integrated lumen that runs distally along the tubing, in which a suction catheter can be inserted to provide access to the stomach for evacuation of the gastric contents.

LarySeal™ Pro is intended to be used to provide oxygen or anaesthetic gases to a patient during surgery and in the pre-hospital setting. The design also allows access to the gastric contents and acts as a conduit for intubation.

LarySeal™ Pro is used to maintain a patent airway in both spontaneously breathing and ventilated patients during anaesthesia and emergency situations. They are a less invasive alternative to an endotracheal tube, particularly in those patients who are at low risk of aspiration, undergoing a short procedure and not given a muscle relaxant.

This device is intended on being used by trained professionals in hospital and emergency environments.

## PRECAUTIONS/WARNINGS

- Do not use if package has been previously opened or damaged, or after the expiry date.
- Keep at ambient storage conditions and avoid exposure to ultraviolet light during storage.
- Before use, visually inspect the device for damage, discolouration or anything that may reduce the functionality of the device.
- Carry out a functional test to ensure there is no occlusion or leaking by inflating the cuff and then deflating the cuff ready for use.
- Lubricate the back of the cuff prior to insertion.
- The position of the LarySeal™ Pro should be reconfirmed after any change in the patient's head or neck position.
- Seat the connector firmly in the adapter on the ventilation equipment to prevent disconnection during use.
- Apply good medical judgement when using this product.
- If the device is incorrectly placed, there is a potential risk of aspiration to the patient.
- In the presence of certain medical gases, the cuff volume or pressure may change. Care should be taken to check the cuff pressure level during the administration of anaesthesia.
- DO NOT apply suction directly to the end of the suction channel.
- DO NOT use in procedures which will involve the use of a laser beam or electro-surgical active electrode in the immediate area or the laryngeal mask airway.
- DO NOT resterilise or reuse.
- Ensure the suction catheter and Endotracheal Tube (ETT) are adequately lubrication prior to insertion.
- DO NOT use non-standard suction catheters or ETT.
- DO NOT use suction catheter or ETT larger than the recommended size applicable to the device.
- DO NOT over inflate the cuff.
- Obstruction may occur if the LarySeal™ Pro becomes dislodged or in incorrectly inserted. Incorrect insertion may cause the epiglottis to block the airway.
- Incorrect placement of the cuff tip into the glottis may mimic bronchospasm.
- The maximum duration of use of the LarySeal™ Pro is 24 hours.

## SPECIFICATION / KEY FEATURES

- Conduit for intubation using an ETT.
- Integrated suction channel for access to gastric contents.
- Mask cuff conforms to the contours of the oropharyngeal area, providing a secure seal for airway management.
- Valve tag prevents pressure build-up in the cuff during transport and storage.
- ISO 5356 compliant 15M connector provides a universal connection.
- ISO 80369-7 compliant valve.
- ISO 17712 compliant supralaryngeal airway device.
- Thin-walled pilot balloon indicates the inflation pressure of the cuff.
- Colour coded pilot balloon indicates mask size.

TABLE 1: LarySeal™ Pro Guidelines on inflation

| Size   | 1      | 1.5    | 2      | 2.5   | 3     | 4     | 5      |
|--|--------|--------|--------|-------|-------|-------|--------|
| PATIENT WEIGHT (kg)                            | <5     | 5-10   | 10-20  | 20-30 | 30-50 | 50-70 | 70-100 |
| MAX CUFF INFLATION VOLUME (ml)                 | 4      | 7      | 10     | 14    | 20    | 30    | 40     |
| MAX CUFF INFLATION PRESSURE (ml)               | 60     | 60     | 60     | 60    | 60    | 60    | 60     |
| MAX ET TUBE                                    | 3.5    | 4      | 5      | 5.5   | 7.5   | 8.0   | 8.0    |
| MAX SUCTION CATHETER (Fr)                      | 8      | 10     | 12     | 14    | 16    | 16    | 16     |
| INTERNAL VOLUME (ml)                           | 6.5    | 7.5    | 14.5   | 20    | 25    | 30    | 30     |
| PRESSURE DROP AT 60L/MIN (cm H <sub>2</sub> O) | <1.0** | <1.0** | <1.0** | <1.0* | <1.0* | <1.0* | <1.0*  |
| MIN INTERDENTAL GAP (mm)                       | 10     | 11     | 13     | 14    | 30    | 30    | 34     |
| NOMINAL LENGTH OF INTERNAL PATHWAY (cm)        | 10     | 11.5   | 11.5   | 13    | 14    | 15.5  | 15.5   |

\*Pressure drop measured at 15L/min

\*\* Pressure drop measured at 30L/min

### PRODUCT RANGE & PART NUMBERS

| Product Code | Product Description                          | Pilot Balloon Colour |
|--------------|--|----------------------|
| 038-94-810   | LarySeal™ Pro Laryngeal Mask Airway Size 1   | White                |
| 038-94-815   | LarySeal™ Pro Laryngeal Mask Airway Size 1.5 | Blue                 |
| 038-94-820   | LarySeal™ Pro Laryngeal Mask Airway Size 2   | Light Green          |
| 038-94-825   | LarySeal™ Pro Laryngeal Mask Airway Size 2.5 | Dark Green           |
| 038-94-830   | LarySeal™ Pro Laryngeal Mask Airway Size 3   | Orange               |
| 038-94-840   | LarySeal™ Pro Laryngeal Mask Airway Size 4   | Red                  |
| 038-94-850   | LarySeal™ Pro Laryngeal Mask Airway Size 5   | Yellow               |

### MATERIALS

| Component   | Material   |
|---|--|
| Connector   | Acrylonitrile Butadiene Styrene (ABS)  |
| Suction Channel   | Polyvinyl Chloride (PVC)   |
| Epiglottic Flap **  | Polyvinyl Chloride (PVC)   |
| Cuff  | Polyvinyl Chloride (PVC)   |
| Shell   | Polyvinyl Chloride (PVC)   |
| Main Tube   | Polyvinyl Chloride (PVC)   |
| Inflation Line  | Polyvinyl Chloride (PVC)   |
| Pilot Balloon   | Polyvinyl Chloride (PVC)   |
| Non return Valve <ul style="list-style-type: none"> <li>• Valve Body</li> <li>• Spring</li> <li>• Valve Rod</li> <li>• Seal Ring</li> </ul> | Polyvinyl Chloride (PVC)<br>Stainless Steel<br>Acrylonitrile Butadiene Styrene (ABS)<br>Silicone |

\*\*Sizes 3-5 only

**LATEX CONTENT**

Flexicare's LarySeal™ Pro does not contain natural rubber latex..

**DEHP CONTENT**

Flexicare's LarySeal™ Pro does not contain phthalate DEHP.

**SINGLE-USE**

Flexicare's LarySeal™ Pro is single use only.

**STERILITY**

Flexicare's LarySeal™ Pro is supplied sterile by Ethylene Oxide (EO) gas.

**STORAGE**

Store in a cool, dry place out of direct sunlight.

**SHELF LIFE**

Flexicare declares a shelf life of 5 years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

**DISPOSAL CONSIDERATIONS**

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

**PACKAGING MATERIALS**

Primary – Polyvinyl Chloride (PVC) and Tyvek

Secondary – Cardboard