smiths medical

portex[®]

Tracheal Tubes

en Instructions for use

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sv Bruksanvisning

da Brugsanvisning

no Bruksanvisning

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Portex® Tracheal Tubes

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® products:

100/105/025-100	Tracheal Tube, Ivory PVC, Oral/Nasal, Uncuffed, sizes 2.5-10.0mm I.D.
100/111/020-110	Tracheal Tube, Siliconised PVC, Oral/Nasal, Uncuffed, sizes 2.0-11.0mm I.D.
100/112/025-050	Tracheal Tube, Paediatric, Oral/Nasal, Uncuffed, Extra long, Siliconised PVC, Croup tube, sizes 2.5-5.0mm I.D.
100/127/025-100	Tracheal Tube, Clear PVC, Murphy Eye, Oral/Nasal, Uncuffed, sizes 2.5-10.0mm I.D.
100/141/025-110	Tracheal Tube, Siliconised PVC, Murphy Eye, Oral/Nasal, Uncuffed, sizes 2.5-11.0mm I.D.
100/166/050-100	Tracheal Tube, Siliconised PVC, Oral/Nasal, Profile® Soft-Seal® (Low Pressure/High Volume) Cuff, sizes 5.0-10.0mm I.D.
100/179/050-090	Tracheal Tube, Ivory PVC, Nasal, Profile® Soft-Seal® (Low Pressure/High Volume) Cuff, sizes 5.0-9.0mm I.D.
100/190/050*	Microlaryngeal Tube, Ivory PVC, Oral/Nasal, Profile® (Low Pressure/High Volume) Cuff, size 5.0mm I.D.
100/199/050-100	Tracheal Tube, Clear PVC, Murphy Eye, Oral/Nasal, Profile® Soft-Seal® (Low Pressure/High Volume) Cuff, sizes 5.0-10.0mm I.D.

*Not for sale in USA

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

DESCRIPTION

A range of sterile, single use tracheal tubes for oral or nasal intubation intended for airway management. Manufactured from clear polyvinyl chloride (PVC), siliconised PVC or Ivory PVC, incorporating the following features:

- Thermosensitive materials with sufficient initial rigidity for intubation which then conforms to the individual patient's respiratory tract at body temperature ensuring minimum trauma.
- Radio-opaque Blue Line® to confirm correct tube placement by X-Ray.

All Portex® tracheal tubes are packed with a 15mm connector conforming to ISO 5356, ISO 7228.

Stricle unless the unit pack is opened or damaged.

INDICATIONS FOR USE

Portex® tracheal tubes are intended for oral and/or nasal intubation for airway management.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- When cutting the tube to length, do not cut close to the joint between the inflation line and the main body of the tracheal tube as blockage or leakage of the inflation system may occur.
- Do not use lubricants to assist the insertion of the 15mm connector into the tracheal tube as disconnection may occur.
- Prior to removal of cuffed tracheal tubes, all air must be completely removed from the cuff to prevent trauma to the airway.
- If the tracheal tube is lubricated prior to insertion, ensure the lubricant does not occlude the lumen of the tube preventing ventilation of the patient. Anaesthetic jellies and lubricants have been associated with complete or partial occlusion of tracheal tubes.
- Do not gauge the adequacy of cuff inflation on the basis of the volume of air introduced, or the level of air resistance felt during inflation.
- During anaesthesia, nitrous oxide may diffuse into the cuff causing an increase in cuff volume.
- Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).
- This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: paediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: www.smiths-medical.com/phthalates

PRECAUTIONS

- The security of all breathing system connectors should be checked when the breathing circuit is established and frequently thereafter. Disconnection may be facilitated with the use of a disconnection wedge (100/555/100) – not available in USA.
- Patients should be adequately humidified to minimise encrustation of the tracheal tube lumen and prevent tracheal mucosal damage.
- The patency of the tracheal tube lumen must be assured by regular suctioning. Check routinely and replace as required to maintain a patent airway. Maximum recommended period of use 30 days.
- Cuff pressure and volume should be monitored and adjusted routinely.
- Devices used in/or during inflation of the cuff must be clean and free from all foreign matter. The inflation device must be removed from the inflation valve immediately after use.
- Guard against cuff damage by avoiding contact with sharp edges. If cuff is damaged, patient should be reintubated and the damaged tube discarded.
- The inflation line valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from the area being scanned.
- Repositioning and movement of the in-situ tracheal tube, while the cuff is inflated, should be avoided.
- In the event that unusual positioning of the head or neck is to be required following intubation, use of a reinforced tube should be considered to avoid the potential for kinking.
- Following intubation, the tracheal tube should be properly secured to help eliminate undesirable movement.
- Follow standard infection control procedures as specified by the Center for Disease Control and Prevention (USA), or local equivalent.

12 The use of topical aerosol anaesthetic agents has been associated with the formation of pinholes in PVC cuffs.

13 Federal (USA) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

(Instructions regarding cuffs and inflation system apply only where these features are present)

- The integrity of the cuff and inflation system should be checked prior to insertion.
- Firmly insert the 15mm connector provided into the tracheal tube.
- Intubate the patient orally or nasally following currently accepted medical techniques. The intubation depth mark may be used as a guide for correct placement of the tube.
- Inflate the cuff with the minimum amount of air to provide an effective seal.

CUFF RESTING DIAMETER (mm)

Tracheal Tube I.D.	5.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0	9.5	10.0	10.5	11.0
Traditional	12	13	15	15	16	17	18	19	21	22	23	24	25
Profile®	15	17	19	20	22	23	25	26	28	30	32	–	–
Profile® Soft-Seal®	17	17	23	23	30	30	30	30	30	30	34	–	–

Blue Line, Profile, Soft-Seal, Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2019 Smiths Medical. All rights reserved.

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en English

Caution

Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 0434A

Do Not Reuse

Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
ISO – 15221 Symbols – Used with Medical Device Labelling, ISO 7000: 1051

Catalogue Number

Indicates the manufacturer's catalogue number so that the medical device can be identified
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 2493

Batch Code

Indicates the manufacturer's batch code so that the batch or lot can be identified
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 2492

Date of Manufacture

Indicates the date when the medical device was manufactured YYYY-MM-DD
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 2497

Use by

Indicates the date after which the medical device is not to be used YYYY-MM-DD
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 2607

Manufacturer

Indicates the medical device manufacturer
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 3082

Not made with natural rubber latex

Large X indicates natural rubber latex (NRL) is not used as a material in the medical product or product container.
ISO-15223-1 Symbols-Used with Medical Device Labelling Annex B.2, ISO 7000: None

Contains or presence of phthalate

Contains or presence of phthalate: bis (2-ethoxyethyl) phthalate (DEHP)
EN 15986:2011 Symbols – Used with Medical Device Labelling – Requirements for labelling of medical devices containing phthalates. Derived from ISO 7000: 2725

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Authorised Representative in the European Community

Indicates the Authorized representative in the European Community.
ISO-15223-1 Symbols - Used with Medical Device Labelling, ISO 7000: None

Sterilized using ethylene oxide

Indicates a medical device that has been sterilized using ethylene oxide
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 2701

Do not use if packaged is damaged

Indicates a medical device that should not be used if the package has been damaged or opened
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 2606

Keep dry

Indicates a medical device that needs to be protected from moisture
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 0626

Keep away from sunlight

Indicates a medical device that needs protection from light sources
ISO – 15223-1 Symbols – Used with Medical Device Labelling, ISO 7000: 0624

Prescription Only

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician
21CFR 801.109 Code Federal Regulations, ISO 7000: None

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Blue Line, Profile, Soft-Seal, Portex, designmerket Portex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet © angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.