









LarySeal™ Pro

LARYNGEAL MASK AIRWAY INSTRUCTIONS FOR USE

EN	Laryngeal Mask Airway	KO	Laryngeal 마스크 기도
PT-BR	Máscara Laringea	LV	Laringālā maska elpceļiem
BG	Ларингеална маска	LT	Gerklų kaukės kvėpavimo takas
HR	Dovod zraka u laringealnu masku	NO	Larynxmaske
CS	Laryngeální maska	PL	Maska krtaniowa
DA	Laryngal luftvejsmaske	PT	Máscara laringea das vias aéreas
NL	Larynxmaskerbeademring	RO	Mască respiratorie laringiană
ET	Kõrimask (LMA)	RU	Ларингеальная маска
FI	Kurkunpäämaskin ilmaputki	SR	Тубус ларингеалне маске
FR	Masque laryngé	ZH	喉罩气道
DE	Larynxmaske für Atemwege	SK	Hrtanová dýchacia cesta pre masku
EL	Αεραγωγός Λαρυγγικής Μάσκας	ES	Máscara laringea
HU	Laringeális maszk	SV	Larynxmask för fri luftväg (LMA)
ID	Sungkup Laring	TR	Laringeal Maske Hava Yolu
IT	Maschera Laringea	UK	Mask Thanh Quản
JA	ラリンジアルマスク気道	VI	Laríngeálnya maska
KK	Airway ларингеалды бетпердеси		

PRODUCT KEY

- 15M Connector
- Main Tube
- Cuff
- Suction Channel
- Shell
- Epiglottic Flap (Sizes 3,4 & 5 only)
- Inflation Line
- Pilot Balloon
- Non-return Valve

Description/Intended Use

LarySeal™ Pro is a supraglottic airway device (SAD) intended to be used to provide oxygen or anaesthetic gases to a patient during surgery and in the pre-hospital setting. The design also allows access to the gastric contents and acts as a conduit for intubation.

Precautions

- Do not use if package has been previously opened or damaged, or after the expiry date.
- Keep at ambient storage conditions and avoid exposure to ultraviolet light during storage.
- Before use, visually inspect the device for damage, discolouration or anything that may reduce the functionality of the device.
- Carry out a functional test to ensure there is no occlusion or leakage by inflating the cuff (refer to Figure 1) and then deflating the cuff ready for use (Figure 2).
- Lubricate the back of the cuff prior to insertion.
- The position of the LarySeal™ Pro should be reconfirmed after any change in the patient's head or neck position.
- Seat the connector firmly in the adapter on the ventilation equipment to prevent disconnection during use.

Warnings & Cautions

- Apply good medical judgment when using this product.
- If the device is incorrectly placed, there is a potential risk of aspiration to the patient.
- In the presence of certain medical gases, the cuff volume or pressure may change. Care should be taken to check the cuff pressure level during the administration of anaesthesia.
- Do not reuse this device. Reuse may result in hazards such as, but not limited to; cross contamination, damage to device and/or suboptimal performance
- Do not apply suction directly to the end of the suction channel.
- Do not use in procedures which will involve the use of a laser beam or electro-surgical active electrode in the immediate area of the laryngeal mask airway.
- Do not resterilise or reuse.
- Ensure the suction catheter and Endotracheal Tube (ETT) are adequately lubricated prior to insertion.
- Incorrect positioning of the LarySeal™ Pro into the patient's airway, may result in ETT not being able to be placed into the trachea
- Do not use non-standard suction catheters or ETT.
- Do not use suction catheter or ETT larger than the recommended size applicable to the device (refer to Table 1).
- Do not over inflate the cuff.
- Obstruction may occur if the LarySeal™ Pro becomes dislodged or is incorrectly inserted. Incorrect insertion may cause the epiglottis to block the airway.
- Incorrect placement of the cuff tip into the glottis may mimic bronchospasm.
- The maximum duration of use of the LarySeal™ Pro is 24 hours.
- Apply good medical judgement when choosing when to remove the LarySeal™ Pro from patient.
- After use, dispose of in accordance with local policies and / or as advised by your local healthcare provider.

Table 1: LarySeal™ Pro Guidelines on Inflation

Size	1	1.5	2	2.5	3	4	5
PATIENT WEIGHT (kg)	<5	5-10	10-20	20-30	30-50	50-79	70-100
MAX CUFF INFLATION VOLUME (ml)	4	7	10	14	20	30	40
MAX CUFF INFLATION PRESSURE (cm H ₂ O)	60	60	60	60	60	60	60
MAX ET TUBE	3.5	4	5	5.5	7.5	8	8
MAX SUCTION CATHETER (Fr)	8	10	12	14	16	16	16
INTERNAL VOLUME (ml)	6.5	7.5	14.5	20	25	30	30
PRESSURE DROP AT 60L/MIN (cm H ₂ O)	<1.0*	<1.0*	<1.0**	<1.0	<1.0	<1.0	<1.0
MIN INTERDENTAL GAP (mm)	10	11	13	14	30	30	34
NOMINAL LENGTH OF INTERNAL PATHWAY (cm)	10	11.5	11.5	13	14	15.5	15.5

*Pressure drop measured at 15L/min

** Pressure drop measured at 30L/min

Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State of the user and/or patient.

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Flexicare Medical Limited