

SPECIFICATIONS
ROMED MEDICAL UNDERPADS,
WITHOUT AND WITH ADHESIVE TAPE
REF: UP60X60, UP60X90, UP60X60-AT, UP60X90-AT



Medical device conform MDR (EU) 2017/745

INTENDED USE

The Romed disposable medical underpads are used to collect moisture and wetness from lying patients in hospitals or at nursing homes, for example postoperatively, with incontinency or after giving birth. The versions with adhesive tape can be fixed to a surface.

MATERIAL

Spun bond nonwoven fabric top layer, moderate absorbent fabric mid layer with superabsorbent polymer (SAP) and waterproof polyethylene (PE) back layer.

Latex free, fragrance free.

DIMENSIONS

Reference no.	Weight	Length	Width	SAP content	Absorbency
UP60X60	50 gram	60 cm	60 cm	4 gram	1000-1100ml
UP60X90	60 gram	60 cm	90 cm	5 gram	1500-1600ml
UP60X60-AT	50 gram	60 cm	60 cm	4 gram	1000-1100ml
UP60X90-AT	60 gram	60 cm	90 cm	5 gram	1500-1600ml

SHELF LIFE

3 years

PACKING

10 pieces in a polybag and 20 polybags (=200 pieces) in an export carton.