



SHEET
WBS – EMERGENCY BLANKET
REF: WBS0003

DATE:
12.2023

RV N°: 3

EMERGENCY BLANKET

1. **DESCRIPTION**

The product WBS EMERGENCY blanket is a disposable self-heating medical device for first aid, which can be used in all emergency situations where it is essential to prevent or treat hypothermia and control body temperature.

2. **INTENDED USE**

The WBS EMERGENCY blanket is a heat source available in a short time that allows you to effectively manage body normothermia and can be used in all cases where it is essential to prevent or treat hypothermia and control body temperature. The device is for single use.

3. **COMPOSITION OF THE PATCHES**

Iron Powder (CAS 7439-89-6) 45~55%
Water (CAS 7732-18-5) 25~31%
Vermiculite (CAS 1318-00-9) 7~8%
Activated Carbon (CAS 64365-11-3) 7~8%
Sodium Chloride (CAS 7647-14-5) 3~5%
Sodium Polyacrylate (CAS 9003-04-7) 1~2%

4. **PACKAGING**

The device is individually vacuum-packed in polyamide and polyethylene sachets.

5. **DISPOSAL**

For recycling, consult the table and check the regulations of your municipality

IFU & BOXES	PAPER	PAP20	
BAG (VACUUM-PACKED)	PLASTIC	PET1	
BLANKET	UNDIFFERENTIATED	PI	
PATCH	UNDIFFERENTIATED	PI	

6. **PRESERVATION**

Store in a cool, dry place, away from heat sources and in tightly closed packages.

ZETA s.r.l.

www.zetagroupsrl.com

Legal Office: Via Loggia dei Pisani 25 -80133- Naples-Italy

Office: Via Atellana 65- 80022 - Arzano - Italy

Tel.: -39 81 64 39 36 **Mail:** info@zetagroupsrl.com

P.IVA: 05712331213

Rea: NA770814

 <p>zeta srl vicini al tuo benessere</p>	<p>SHEET WBS – EMERGENCY BLANKET REF: WBS0003</p>	<p>DATE: 12.2023</p>	<p>RV N°: 3</p>
--	---	--------------------------	-----------------

7. CLASSIFICATION

Class I non-sterile disposable medical device CE marked according to EU Regulation 2017/745.

8. CND

T030399 - HEATING COVERS – OTHER


9. RDM

2426717

10. BASIC UDI-DI

803354931WBSL5

11. SPECIFICATIONS

<u>EMERGENCY BLANKET</u>	<u>Dimension</u>	<u>Composition:</u>	<u>How to use:</u>	<u>Destination:</u>
	<p><u>Open:</u> 150*90 cm</p> <p><u>Lock:</u> 27*31 cm</p>	<p>Two layers of SMS polypropylene non-woven fabric with pockets that contain 12 self-heating patches.</p>	<p>Remove the device from the packaging, wait 15 minutes, after exposure to air, the Heating Patches will start to generate heat.</p>	<p>All emergency and first aid situations where hypothermia needs to be checked</p>
	<p><u>Weight:</u> 648.4gr</p>	<p><u>Duration:</u> 14h+</p>		<p><u>UDI Code OF:</u> 08033549310530</p>
		<p><u>Interval Temperature:</u> 48°C<<T<<52°C</p>		<p><u>REF:</u> WBS0003</p>

ZETA s.r.l.

www.zetagroupsrl.com

Legal Office: Via Loggia dei Pisani 25 -80133- Naples-Italy

Office: Via Atellana 65- 80022 - Arzano - Italy

Tel.: -39 81 64 39 36 **Mail:** info@zetagroupsrl.com

P.IVA: 05712331213

Rea: NA770814



SHEET
WBS – EMERGENCY BLANKET
REF: WBS0003

DATE:
12.2023

RV N°: 3

ZETA s.r.l.

www.zetagroupsrl.com

Legal Office: Via Loggia dei Pisani 25 -80133- Naples-Italy

Office: Via Atellana 65- 80022 - Arzano - Italy

Tel.: -39 81 64 39 36 **Mail:** info@zetagroupsrl.com

P.IVA: 05712331213

Rea: NA770814